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CONTENTS

1. AREA OF APPLICATION
2. FEATURES
3. SAFETY
4. INSTALLATION AND CONNECTING TO
CYLINDER / BUNDLE
5. START-UP AND NORMAL USE
6. MAINTENANCE, CLEANING
AND RECYCLING

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IMPORTANT

In order to maintain the quality of our product in the best safety conditions throughout its period of use, we strongly recommend you to read this brochure and the NT0100 attached to this product.

It is imperative to follow the advices they contain, and that applies to every paragraph of this manual.

Ind.	Modifications	Nom	Date	Ind.	Modifications	Name	Date
H	Edition modifiée	CCH	01/21	H	Modified issue	CCH	01/21
I	Edition modifiée	MOG	05/22	I	Modified issue	MOG	05/22
J	Edition modifiée :	DDA	11/23	J	Modified issue:	DDA	11/23

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§ 1: Addition of the high-pressure flexible fittings standard
§ 2: General update

52, rue de Lagny
93100 MONTREUIL - FRANCE
Tél : +33 (0)1.41.72.90.00
Fax : +33 (0)1.41.72.90.01
E-mail : cahouet@cahouet.com
www.cahouet.com

NT 0572 J
FR/EN

FLEXIBLES HAUTE PRESSION HIGH PRESSURE FLEXIBLE HOSES

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1. AREA OF APPLICATION

High pressure flexible hoses are designed to be used with all medical gases, included oxygen. Our flexible hoses are part of a medical gas pipeline system in accordance with the standard EN ISO 7396-1 and meet the requirements of the high-pressure flexible hose fittings standard NF EN ISO 21969 for use with the medical gas systems.

2. FEATURES

- Nominal pressure (P1): 200 bar
- Testing pressure: > 375 bar
- Flow coefficient :
 - cylinder hose: Kv = 3,97
 - bundle hose: Kv = 8,75
- Flexible tube:
 - Inner tube Ø: 6 or 10 mm
 - Internal tube: corrugated stainless steel
 - Reinforcement: stainless steel braids
- Safety cable: stainless steel
- Operating or storage temperature: - 20°C to + 60°C
- Length:
 - cylinder hose: 0.5 to 3 m
 - bundle hose: 0.5 to 4 m
 - flexible hose with yoke connection: 0.5 to 1.5 m
- Inlet connection: according to gas (standard/country)
- Outlet connection: SI 21,7x1,814 RH, M 20x150 RH or other connections (standard/country) with an opening for sealing or gas specificity.
- Free of grease for oxygen use.

3. SAFETY

This equipment works with gas under pressure. You must read and follow our general rules for installation and operation supplied with the equipment (NT0100).

- Never let a flexible hose hanging down from one of its ends. Penetrations of particles while manipulating are responsible for most breakdowns.
- Non connected ends of HP flexible hoses must be protected by a plug.
- The flexible hose should never touch the floor or any other element.
- The installation date must absolutely be marked on the hose plate by the installer.
- In case of non gas specific flexible hoses, seal the connection hose/collector.

4. INSTALLATION AND CONNECTING TO CYLINDER / BUNDLE

- Check the presence and the good condition of the sealing components. They must be without damage or scratches in order to avoid any risk of leakage.
- Both ends of the flexible hose have to be in the same vertical plane.
- The radius of curvature of the hose must be higher than 110 mm (Ø 6) and 150 mm (Ø 10).

Examples of installation:

The tightening torque for JIC connections must conform the ISO 8434-2 standard (ex. 15 N.m maxi pour 7/16").



NT 0572 J
FR/EN


FLEXIBLES HAUTE PRESSION HIGH PRESSURE FLEXIBLE HOSES

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- Remove the protection caps and check that the seals are fitted at both ends of the hose except in case of metal/metal tightness.
- In case of an O-ring connection:
 - 1) Check that the O-ring is well positioned in its groove to avoid any cut.
 - 2) The sealing O-ring must be replaced each time disconnection.
- At the collector side, use two wrenches to prevent the non return valves or the intermediate fitting from turning.
- When tightening the flexible hoses on to the collector's inlet ports, the torque should be limited to 50 N.m for metal/metal tightness and about 8 to 10 N.m for O-ring tightness types.
- Fasten the safety cables through the holes provided on the panel of the change-over and/or around the collectors and/or onto a firm anchor point.
- Branch and handtighten on the cylinders/bundles valves the high-pressure hoses with care (around 8 to 10 N.m).
- Tie up the safety cable loop around the cylinder / bundle valve.
- Seal the connection collector/hose to avoid any unexpected intervention.

5. START-UP AND NORMAL USE

For testing, start-up and conformity of the installation, refer to the standard ISO 7396-1.

- Close the upstream and downstream pipeline system.
- Connect both ends of the flexible hose.
- Slowly open and in 2 steps the gas bundle valve or the cylinder valve or the upstream pipeline system :
 - 1) Slowly and partially open.
Wait for a pressure stability.
 - 2) Then slowly finish the opening of the valve and turn it 1/4 turn backwards.
- Check the tightness with a non-fat foam reagent.

Purge the hose and start by repeating the above-mentioned procedure.

ATTENTION: Drain the flexible hose before any disconnecting and before any replacement with a hose of identical features.

NT 0572 J
FR/EN

FLEXIBLES HAUTE PRESSION HIGH PRESSURE FLEXIBLE HOSES

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6. MAINTENANCE, CLEANING AND RECYCLING

The safe use of flexible hoses requires the following regular precautions:

- Control of the shaft: flexible hoses with more than 2 filaments of the braid broken per cm^2 must be changed.
- Control the condition of the threads and check the lack of corrosion during each cylinder/bundle change.
- Change the sealing ring after each dismounting.

Maintenance kits (O-Rings or washers) are available upon request indicating the references of the concerned hoses.

- A HP flexible hose must be discarded after 5 years from the commissioning date or 2000 cycles of connection/disconnection carried out by competent and trained personnel. The lower value shall be retained.
- Change the hose if the initial flow is no longer attained.
- When changing a HP flexible hose, only use guaranteed original HP flexible hoses, which are the only ones compatible with your installation.

NT 0572 J
FR/EN